


SCHOTT LABORATORY GLASSWARE PRODUCT LISTING


Beaker, low form, with spout, DURAN borosilicate

<i>Cat. No</i>	<i>capacity ml</i>	<i>diameter d, mm</i>	<i>height h, mm</i>	<i>pack size</i>
21 106 07	5	22	30	10
21 106 08	10	26	35	10
21 106 14	25	34	50	10
21 106 17	50	42	60	10
21 106 24	100	50	70	10
21 106 29	150	60	80	10
21 106 36	250	70	95	10
21 106 41	400	80	110	10
21 106 48	600	90	125	10
21 106 53	800	100	135	10
21 106 54	1000	105	145	10
21 106 63	2000	132	185	10
21 106 68	3000	152	210	4
21 106 73	5000	170	270	4
21 106 86	10000	217	350	1




Beaker, tall form, with spout, DURAN borosilicate

<i>Cat. No</i>	<i>capacity ml</i>	<i>diameter d, mm</i>	<i>height h, mm</i>	<i>pack size</i>
21 116 17	50	38	70	10
21 116 24	100	48	80	10
21 116 29	150	54	95	10
21 116 36	250	60	120	10
21 116 41	400	70	130	10
21 116 48	600	80	150	10
21 116 53	800	90	175	10
21 116 54	1000	95	180	10
21 116 63	2000	120	240	10
21 116 68	3000	135	280	2




Beaker, tall form, without spout, DURAN borosilicate

<i>Cat. No</i>	<i>capacity ml</i>	<i>diameter d, mm</i>	<i>height h, mm</i>	<i>pack size</i>
21 117 17	50	38	70	10
21 117 24	100	48	80	10
21 117 29	150	54	95	10
21 117 36	250	60	120	10
21 117 41	400	70	130	10
21 117 48	600	80	150	10
21 117 54	1000	95	180	10



Beaker, heavy wall, with spout, DURAN borosilicate


<i>Cat. No</i>	<i>capacity ml</i>	<i>diameter d, mm</i>	<i>height h, mm</i>	<i>pack size</i>
21 131 24	100	52	85	10
21 131 29	150	54	93	10
21 131 36	250	70	94	10
21 131 44	500	89	124	10
21 131 54	1000	105	160	10



SCHOTT LABORATORY GLASSWARE PRODUCT LISTING


Beaker, heavy wall, with spout, DURAN borosilicate, continued

<i>Cat. No</i>	<i>capacity ml</i>	<i>diameter d, mm</i>	<i>height h, mm</i>	<i>pack size</i>
21 131 63	2000	135	195	10
21 131 68	3000	157	205	4
21 131 73	5000	182	256	1
21 131 86	10000	225	340	1
21 131 88	15000	260	390	1
21 131 91	20000	285	430	1




Beaker, Philips, with spout, DURAN borosilicate

<i>Cat. No</i>	<i>capacity ml</i>	<i>diameter d, mm</i>	<i>height h, mm</i>	<i>pack size</i>
21 141 29	150	59	87	10
21 141 36	250	68	105	10
21 141 44	500	86	142	10




Beaker, Berzelius, DURAN borosilicate

<i>Cat. No</i>	<i>capacity ml</i>	<i>diameter d, mm</i>	<i>height h, mm</i>	<i>pack size</i>
21 126 01	100	50	78	10




Vessel, Bloom test, DURAN borosilicate

<i>Cat. No</i>	<i>diameter d, mm</i>	<i>height h, mm</i>	<i>pack size</i>
21 125 01	59	85	10



Erlenmeyer flask, narrow neck, DURAN borosilicate


<i>Cat. No</i>	<i>capacity ml</i>	<i>diameter d, mm</i>	<i>diameter d1, mm</i>	<i>height h, mm</i>	<i>pack size</i>
21 216 14	25	42	22	75	10
21 216 17	50	51	22	90	10
21 216 24	100	64	22	105	10
21 216 32	200	79	34	131	10
21 216 36	250	85	34	145	10
21 216 39	300	87	34	156	10
21 216 44	500	105	34	180	10
21 216 54	1000	131	42	220	10
21 216 63	2000	166	50	280	2
21 216 68	3000	181	51	310	2
21 216 73	5000	220	51	365	1



SCHOTT LABORATORY GLASSWARE PRODUCT LISTING


Erlenmeyer flask, wide neck, DURAN borosilicate

Cat. No	capacity ml	diameter d, mm	diameter d 1, mm	height h, mm	pack size
21 226 14	25	43	31	70	10
21 226 17	50	51	34	85	10
21 226 24	100	64	34	105	10
21 226 32	200	79	50	131	10
21 226 36	250	85	50	140	10
21 226 39	300	87	50	156	10
21 226 44	500	105	50	175	10
21 226 54	1000	131	50	220	10
21 226 63	2000	153	72	276	10
21 227 68	3000	190	106	285	1
21 227 73	5000	220	108	322	1
21 227 86	10000	285	147	420	1




Erlenmeyer flask with screw-threaded cap, DURAN borosilicate

Cat. No	capacity ml	DIN thread GL	diameter d, mm	height h, mm	pack size
21 803 245	100	25	64	109	10
21 803 365	250	32	85	149	10
21 803 445	500	32	105	180	10
21 803 545	1000	32	131	225	10




Flask with standard ground joint and glass stopper, (iodine flask), DURAN borosilicate

Cat. No	capacity ml	diameter d, mm	height h, mm	pack size
24 192 27	100	64	120	10
24 192 37	250	85	160	10
24 192 46	500	105	195	10
24 192 56	1000	131	235	10



Flask, round bottom, narrow neck, DURAN borosilicate


Cat. No	capacity ml	diameter d, mm	diameter d 1, mm	height h, mm	pack size
21 721 17	50	51	26	95	10
21 721 24	100	64	26	110	10
21 721 36	250	85	34	144	10
21 721 44	500	105	34	175	10
21 721 54	1000	134	42	200	10
21 721 64	2000	166	42	260	10
21 721 68	3000	185	50	260	1
21 721 71	4000	207	50	290	1
21 721 73	5000	223	50	305	1
21 721 77	6000	236	51	355	1
21 721 86	10000	279	65	420	1
21 721 87	12000	295	71	380	1
21 721 91	20000	345	76	515	1



SCHOTT LABORATORY GLASSWARE PRODUCT LISTING


Flask, round bottom, wide neck, DURAN borosilicate

Cat. No	capacity ml	diameter d, mm	diameter d 1, mm	height h, mm	pack size
21 741 17	50	51	34	105	10
21 741 24	100	64	35	110	10
21 741 36	250	85	51	145	10
21 741 44	500	105	50	168	10
21 741 54	1000	131	50	210	10
21 741 55	1000	131	65	210	10
21 741 63	2000	166	76	260	10
21 741 64	2000	166	50	260	10
21 741 68	3000	185	65	260	1
21 741 71	4000	206	76	315	1
21 741 73	5000	223	65	310	1
21 741 76	6000	236	89	355	1
21 741 77	6000	236	65	355	1
21 741 86	10000	279	89	420	1
21 741 91	20000	345	89	520	1




Flask, flat bottom, narrow neck, DURAN borosilicate

Cat. No	capacity ml	diameter d, mm	diameter d 1, mm	height h, mm	pack size
21 711 17	50	51	34	90	10
21 711 24	100	64	22	110	10
21 711 36	250	85	35	140	10
21 711 44	500	105	34	171	10
21 711 54	1000	131	42	190	10
21 711 64	2000	166	42	250	10
21 711 68	3000	185	50	250	1
21 711 71	4000	207	50	300	1
21 711 73	5000	223	50	290	1
21 711 76	6000	237	65	340	1
21 711 86	10000	280	65	400	1




Flask, flat bottom, wide neck, DURAN borosilicate

Cat. No	capacity ml	diameter d, mm	diameter d 1, mm	height h, mm	pack size
21 731 17	50	51	35	100	10
21 731 24	100	64	34	110	10
21 731 36	250	85	50	140	10
21 731 44	500	103	50	170	10
21 731 54	1000	131	50	200	10
21 731 63	2000	166	76	250	10
21 731 64	2000	166	50	250	10



Flask, distillation, DURAN borosilicate


Cat. No	capacity ml	diameter d, mm	height h, mm	side arm d 1 / l mm	pack size
21 243 24	100	63	150	8/200	10
21 243 36	250	85	200	8/200	10
21 243 44	500	105	250	8/200	10
21 243 54	1000	131	300	12/250	10



SCHOTT LABORATORY GLASSWARE PRODUCT LISTING


Flask, distillation, Engler, DURAN borosilicate

<i>Cat. No</i>	<i>capacity ml</i>	<i>diameter d / d 1 mm</i>	<i>height h, mm</i>	<i>side arm d 2 / l mm</i>	<i>pack size</i>
21 653 24	100	66/16	215	6/100	10
21 653 28	125	68/18	215	7/100	10
21 653 29	150	73/16	223	6/100	10




Flask, distillation, Engler, with standard 19/26 ground joint, DURAN borosilicate

<i>Cat. No</i>	<i>capacity ml</i>	<i>diameter d, mm</i>	<i>height h, mm</i>	<i>side arm d 1 / l mm</i>	<i>pack size</i>
24 653 28	125	68	215	7/100	10




Flask, Kjeldahl, DURAN borosilicate

<i>Cat. No</i>	<i>capacity ml</i>	<i>diameter d, mm</i>	<i>height h, mm</i>	<i>neck dia. d 1, mm</i>	<i>pack size</i>
21 231 17	50	51	200	22	10
21 231 24	100	60	200	22	10
21 231 36	250	81	270	34	10
21 231 44	500	101	300	34	10
21 231 51	750	115	340	34	10
21 231 54	1000	126	350	34	10



Glass ceramic laboratory protection plate

<i>Cat. No</i>	<i>dimensions a x b mm</i>	<i>pack size</i>
23 821 53	135 x 135	10
23 821 57	155 x 155	10
23 821 58	175 x 175	10



Quadrupod for glass ceramic laboratory protection plate, chrome-nickel steel

<i>Cat. No</i>	<i>dimensions a x b mm</i>	<i>height h, mm</i>	<i>pack size</i>
29 077 53	135 x 135	220	5
29 077 57	155 x 155	220	5
29 077 58	175 x 175	220	5





Plate holder for glass ceramic laboratory protection plate, chrome-nickel steel


<i>Cat. No</i>	<i>dimensions a x b mm</i>	<i>pack size</i>
29 078 53	135 x 135	5
29 078 57	155 x 155	5
29 078 58	175 x 175	5



SCHOTT LABORATORY GLASSWARE PRODUCT LISTING

Spirit lamp, soda lime, without filler tubulature, with ground-on cap, without socket and wick

<i>Cat. No</i>	<i>capacity ml</i>	<i>diameter d, mm</i>	<i>height h, mm</i>	<i>pack size</i>
23 400 24	100	75	103	10



Spirit lamp, soda lime, without filler tubulature, with ground-on cap, with socket and wick

<i>Cat. No</i>	<i>capacity ml</i>	<i>diameter d, mm</i>	<i>height h, mm</i>	<i>pack size</i>
23 400 245	100	75	103	10

Spirit lamp, soda lime, with filler tubulature and stopper, with ground-on cap, without socket and wick

<i>Cat. No</i>	<i>capacity ml</i>	<i>diameter d, mm</i>	<i>height h, mm</i>	<i>pack size</i>
23 401 24	100	75	103	10



Spirit lamp, soda lime, with filler tubulature and stopper, with ground-on cap, with socket and wick

<i>Cat. No</i>	<i>capacity ml</i>	<i>diameter d, mm</i>	<i>height h, mm</i>	<i>pack size</i>
23 401 245	100	75	103	10

Wick for spirit lamp

<i>Cat. No.</i>	<i>pack size</i>
29 402 00	1




Socket for spirit lamp

<i>Cat. No.</i>	<i>pack size</i>
29 403 00	1



Evaporating dish, DURAN borosilicate


<i>Cat. No</i>	<i>capacity ml</i>	<i>diameter d, mm</i>	<i>height h, mm</i>	<i>pack size</i>
21 301 32	15	50	25	10
21 301 34	45	60	30	10
21 301 38	60	70	35	10
21 301 41	90	80	45	10
21 301 44	170	95	55	10
21 301 49	320	115	65	10
21 301 54	600	140	80	10
21 301 59	1500	190	100	1
21 301 63	2500	230	130	1



SCHOTT LABORATORY GLASSWARE PRODUCT LISTING


Crystallising dish, DURAN borosilicate

<i>Cat. No</i>	<i>capacity ml</i>	<i>diameter d, mm</i>	<i>height h, mm</i>	<i>pack size</i>
with spout				
21 311 24	20	40	25	10
21 311 32	40	50	30	10
21 311 34	60	60	35	10
21 311 38	100	70	40	10
21 311 41	150	80	45	10
21 311 44	300	95	55	10
21 311 49	500	115	65	10
21 311 54	900	140	75	10
21 311 59	2000	190	90	10
21 311 63	3500	230	100	1




Crystallising dish, DURAN borosilicate

<i>Cat. No</i>	<i>capacity ml</i>	<i>diameter d, mm</i>	<i>height h, mm</i>	<i>pack size</i>
without spout				
21 313 24	20	40	25	10
31 313 32	40	50	30	10
21 313 34	60	60	35	10
21 313 38	100	70	40	10
21 313 41	150	80	45	10
21 313 44	300	95	55	10
21 313 49	500	115	65	10
21 313 54	900	140	75	10
21 313 59	2000	190	90	10
21 313 63	3500	230	100	1




Weighing bottle, with grip stopper, DURAN borosilicate

<i>Cat. No</i>	<i>capacity ml</i>	<i>diameter d, mm</i>	<i>height h, mm</i>	<i>pack size</i>
low form				
24 210 13	5	28	25	10
24 210 23	15	38	30	10
24 210 32	30	54	30	10
24 210 41	80	85	30	10
tall form				
24 211 13	10	28	40	10
24 211 18	20	32	50	10
24 211 23	45	38	70	10
24 211 24	70	44	80	10



Watch glass, with fused edge

<i>Cat. No.</i>	<i>d mm</i>	<i>pack size</i>
DURAN borosilicate		
21 321 24	40	10
21 321 32	50	10
21 321 34	60	10
21 321 41	80	10
21 321 46	100	10



SCHOTT LABORATORY GLASSWARE PRODUCT LISTING

Watch glass, with fused edge

<i>Cat. No.</i>	<i>d</i> <i>mm</i>	<i>pack size</i>
-----------------	-----------------------	------------------


DURAN borosilicate, continued

21 321 52	125	10
21 321 57	150	10
21 321 61	200	10
21 321 66	250	1



soda-lime


23 321 24	40	10
23 321 32	50	10
23 321 34	60	10
23 321 38	70	10
23 321 41	80	10
23 321 43	90	10
23 321 46	100	10
23 321 51	120	10
23 321 52	125	10
23 321 57	150	10
23 321 61	200	10
23 321 66	250	1



Cylinder, with knobbed lid, polished rim, DURAN borosilicate

<i>Cat. No.</i>	<i>d</i> <i>mm</i>	<i>h</i> <i>mm</i>	<i>pack size</i>
-----------------	-----------------------	-----------------------	------------------


24 205 01	80	80	10
24 205 03	100	100	10
24 205 05	120	120	1
24 205 10	150	150	1
24 205 21	210	210	1
24 205 32	260	260	1



Cylinder, multi-purpose, round base, rough ground rim, DURAN borosilicate

<i>Cat. No</i>	<i>capacity</i> <i>ml</i>	<i>diameter</i> <i>d, mm</i>	<i>height</i> <i>h, mm</i>	<i>pack size</i>
----------------	------------------------------	---------------------------------	-------------------------------	------------------


21 398 21	300	50	150	10
21 398 34	250	40	200	10
21 398 36	570	60	200	10
21 398 46	700	60	250	10
21 398 52	380	40	300	10
21 398 53	600	50	300	10
21 398 68	500	40	400	10
21 398 74	2000	80	400	10
21 398 77	1500	65	450	10
21 398 80	1000	50	500	1



Bell jar, with knob, DURAN borosilicate

<i>Cat. No</i>	<i>height</i> <i>h 1, mm</i>	<i>height</i> <i>h, mm</i>	<i>diameter</i> <i>d, mm</i>	<i>pack size</i>
----------------	---------------------------------	-------------------------------	---------------------------------	------------------

24 460 59	50	250	185	1
24 460 66	50	255	260	1
24 460 69	50	300	315	1



SCHOTT LABORATORY GLASSWARE PRODUCT LISTING


Bell jar, with NS 34/35 neck bore, DURAN borosilicate

Cat. No	height h 1, mm	height h, mm	diameter d, mm	pack size
24 465 59	50	250	185	1
24 465 61	50	300	215	1
24 465 69	50	500	315	1




Sedimentation cone, Imhoff, graduated, DURAN borosilicate

Cat. No	capacity ml	max. dia. d, mm	height h, mm	pack size
21 401 54	1000	120	465	1




Sedimentation cone, Imhoff, graduated, with stopcock, DURAN borosilicate

Cat. No	capacity ml	max. dia. d, mm	height h, mm	pack size
21 402 454	1000	120	465	1




Sedimentation cylinder, graduated, with NS 34/35 stopper, DURAN borosilicate

Cat. No	capacity ml	diameter d, mm	height h, mm	pack size
21 619 54	500	54	395	1




Laboratory bottle, with DIN thread, without cap or pouring ring, DURAN borosilicate

Cat. No	capacity ml	DIN thread GL	diameter d, mm	height h, mm	pack size
21 801 14	25	25	36	70	10
21 801 17	50	32	46	87	10
21 801 24	100	45	56	100	10
21 801 36	250	45	70	138	10
21 801 44	500	45	86	176	10
21 801 54	1000	45	101	225	10
21 801 63	2000	45	136	260	10
21 801 73	5000	45	183	330	10
21 801 86	10000	45	227	410	1
21 801 88	15000	45	268	445	1
21 801 91	20000	45	288	505	1



Laboratory bottle, DIN thread, with PP cap & PP pouring ring, DURAN borosilicate

Cat. No	capacity ml	DIN thread GL	diameter d, mm	height h, mm	pack size
21 801 145	25	25	36	74	10
21 801 175	50	32	46	91	10
21 801 245	100	45	56	105	10
21 801 365	250	45	70	143	10
21 801 445	500	45	86	181	10
21 801 545	1000	45	101	230	10
21 801 635	2000	45	136	265	10
21 801 735	5000	45	182	335	2




SCHOTT LABORATORY GLASSWARE PRODUCT LISTING

Laboratory bottle, DIN thread, with PP cap & PP pouring ring, DURAN borosilicate cont.

<i>Cat. No</i>	<i>capacity ml</i>	<i>DIN thread GL</i>	<i>diameter d, mm</i>	<i>height h, mm</i>	<i>pack size</i>
21 801 865	10000	45	227	415	1
21 801 885	15000	45	268	450	1
21 801 915	20000	45	288	510	1


Laboratory bottle, with DIN thread, without cap or pouring ring, AMBER DURAN borosilicate

<i>Cat. No</i>	<i>capacity ml</i>	<i>DIN thread GL</i>	<i>diameter d, mm</i>	<i>height h, mm</i>	<i>pack size</i>
21 806 14	25	25	36	70	10
21 806 17	50	32	46	87	10
21 806 24	100	45	56	100	10
21 806 36	250	45	70	138	10
21 806 44	500	45	86	176	10
21 806 54	1000	45	101	225	10
21 806 63	2000	45	136	260	10
21 806 73	5000	45	182	330	1
21 806 86	10000	45	227	410	1
21 806 88	15000	45	268	445	1
21 806 91	20000	45	288	505	1




Laboratory bottle, with DIN thread, plastic coated, without cap or pouring ring, DURAN borosilicate

<i>Cat. No</i>	<i>capacity ml</i>	<i>DIN thread GL</i>	<i>diameter d, mm</i>	<i>height h, mm</i>	<i>pack size</i>
21 805 24	100	45	56	100	10
21 805 36	250	45	70	138	10
21 805 44	500	45	86	176	10
21 805 54	1000	45	101	225	10
21 805 63	2000	45	136	260	10
21 805 73	5000	45	182	330	1



Laboratory bottle, with DIN thread, pressure resistant to -1/+1.5 bar, without cap or pouring ring, DURAN borosilicate

<i>Cat. No</i>	<i>capacity ml</i>	<i>DIN thread GL</i>	<i>diameter d, mm</i>	<i>height h, mm</i>	<i>pack size</i>
clear glass 21 810 54	1000	45	101	225	10
amber glass 21 816 54	1000	45	101	225	10
plastic coated 21 815 54	1000	45	101	225	



Laboratory bottle, with DIN thread, square, with PP cap and PP pouring ring, DURAN borosilicate


<i>Cat. No</i>	<i>capacity ml</i>	<i>DIN thread GL</i>	<i>diameter d, mm</i>	<i>height h, mm</i>	<i>pack size</i>
21 820 245	100	32	50	109	10
21 820 365	250	45	64	143	10
21 820 445	500	45	78	181	10
21 820 545	1000	45	94	222	10



SCHOTT LABORATORY GLASSWARE PRODUCT LISTING


Screw caps, polypropylene, autoclavable to 140°C steam

<i>Cat. No</i>	<i>colour</i>	<i>DIN thread GL</i>	<i>diameter d, mm</i>	<i>height h, mm</i>	<i>pack size</i>
29 239 13	blue	25	33	19	10
29 239 19	blue	32	40	24	10
29 239 28	blue	45	54	26	10
29 338 280	yellow	45	54	26	10
29 338 286	green	45	54	26	10
29 338 288	grey	45	54	26	10




Pouring ring, polypropylene, autoclavable to 140°C steam

<i>Cat. No</i>	<i>colour</i>	<i>DIN thread GL</i>	<i>height h, mm</i>	<i>pack size</i>
29 242 19	blue	32	4	10
29 242 28	blue	45	4	10
10 899 11	green	45	4	10
10 899 14	grey	45	4	10
10 899 17	yellow	45	4	10




Screw caps, polypropylene, tamper-proof, autoclavable to 140°C steam

<i>Cat. No</i>	<i>colour</i>	<i>DIN thread GL</i>	<i>diameter d, mm</i>	<i>height h, mm</i>	<i>pack size</i>
10 175 26	blue	45	54	26	10




Screw caps, TpCh260, premium, useable range -196°C to + 260°C

<i>Cat. No</i>	<i>colour</i>	<i>DIN thread GL</i>	<i>diameter d, mm</i>	<i>height h, mm</i>	<i>pack size</i>
10 886 79	no colour	45	52	27	5




Pouring ring, TpCh260, usable range -196°C to + 260°C

<i>Cat. No</i>	<i>colour</i>	<i>DIN thread GL</i>	<i>height h, mm</i>	<i>pack size</i>
10 886 78	no colour	45	4	5




Screw cap, PP, with aperture and ptfе membrane, autoclavable to 140°C

<i>Cat. No</i>	<i>colour</i>	<i>DIN thread GL</i>	<i>diameter d, mm</i>	<i>height h, mm</i>	<i>pack size</i>
10 886 55	blue	45	54	28	5



Screw cap, PBT, red, with ptfе protected seal, will withstand hot air sterilisation to 180°C


<i>Cat. No</i>	<i>DIN thread GL</i>	<i>diameter d, mm</i>	<i>height h, mm</i>	<i>pack size</i>
29 240 13	25	33	23	10
29 240 19	32	40	26	10
29 240 28	45	54	28	10



SCHOTT LABORATORY GLASSWARE PRODUCT LISTING


Pouring ring, EFTE, red, will withstand hot air sterilisation to 180°C

Cat. No	DIN thread GL	height h, mm	pack size
29 244 19	32	4	10
29 244 28	45	4	10




Bottle, reagent, wide neck, with standard ground joint, DURAN borosilicate

Cat. No	capacity ml	dia. d, mm	height h, mm	neck	pack size
clear glass					
21 184 17	50	44	79	24/20	10
21 184 24	100	52	97	29/22	10
21 184 36	250	70	132	34/35	10
21 184 44	500	86	163	45/40	10
21 184 54	1000	107	201	60/46	10
21 184 63	2000	133	247	60/46	10
21 184 73	5000	182	358	85/55	1




Bottle, reagent, wide neck, with standard ground joint, DURAN borosilicate

Cat. No	capacity ml	dia. d, mm	height h, mm	neck	pack size
clear glass, continued					
21 184 86	10000	229	443	85/55	1
21 184 91	20000	290	570	85/55	1
amber glass					
21 187 17	50	44	79	24/20	10
21 187 24	100	52	87	29/22	10
21 187 36	250	70	132	34/35	10
21 187 44	500	86	163	45/40	10
21 187 54	1000	107	201	60/46	10
21 187 63	2000	133	247	60/46	10
21 187 73	5000	182	358	85/55	1
21 187 86	10000	229	443	85/55	1
21 187 91	20000	290	570	85/55	1



Bottle, reagent, wide neck, with glass flat-head stopper, DURAN borosilicate


Cat. No	capacity ml	dia. d, mm	height h, mm	neck	pack size
clear glass					
21 185 17	50	44	79	24/20	10
21 185 24	100	52	97	29/22	10
21 185 36	250	70	132	34/35	10
21 185 44	500	86	163	45/40	10
21 185 54	1000	107	201	60/46	10
21 185 63	2000	133	247	60/46	10
21 185 73	5000	182	358	85/55	1
21 185 86	10000	229	443	85/55	1
21 185 91	20000	290	570	85/55	1



SCHOTT LABORATORY GLASSWARE PRODUCT LISTING


Bottle, reagent, wide neck, with glass flat-head stopper, DURAN borosilicate, continued

<i>Cat. No</i>	<i>capacity ml</i>	<i>dia. d, mm</i>	<i>height h, mm</i>	<i>neck</i>	<i>pack size</i>
amber glass					
21 188 17	50	44	79	24/20	10
21 188 24	100	52	97	29/22	10
21 188 36	250	70	132	34/35	10
21 188 44	500	86	163	45/40	10
21 188 54	1000	107	201	60/46	10
21 188 63	2000	133	247	60/46	10
21 188 73	5000	182	358	85/55	1
21 188 86	10000	229	443	85/55	1
21 188 91	20000	290	570	85/55	1




Bottle, reagent, wide neck, with standard ground joint, SODA-LIME

<i>Cat. No</i>	<i>capacity ml</i>	<i>dia. d, mm</i>	<i>height h, mm</i>	<i>neck</i>	<i>pack size</i>
clear glass					
23 184 17	50	44	78	24/20	10
23 184 24	100	54	100	29/22	10
23 184 36	250	72	130	34/24	10
23 184 44	500	89	165	45/27	10
23 184 54	1000	110	200	60/46	10
amber glass					
23 187 17	50	44	78	24/20	10
23 187 24	100	54	100	29/22	10
23 187 36	250	72	130	34/24	10
23 187 44	500	89	165	45/27	10
23 187 54	1000	110	200	60/46	10



Bottle, reagent, wide neck, with glass flat-head stopper, SODA-LIME

<i>Cat. No</i>	<i>capacity ml</i>	<i>dia. d, mm</i>	<i>height h, mm</i>	<i>neck</i>	<i>pack size</i>
clear glass					
23 185 17	50	44	78	24/20	10
23 185 24	100	54	100	29/22	10
23 185 36	250	72	130	34/24	10
23 185 44	500	89	165	45/27	10
23 185 54	1000	110	200	60/46	10
amber glass					
23 188 17	50	44	78	24/20	10
23 188 24	100	54	100	29/22	10
23 188 36	250	72	130	34/24	10
23 188 44	500	89	165	45/27	10
23 188 54	1000	110	200	60/46	10



SCHOTT LABORATORY GLASSWARE PRODUCT LISTING


Bottle, reagent, wide neck, with polyethylene stopper, SODA-LIME

<i>Cat. No</i>	<i>capacity ml</i>	<i>dia. d, mm</i>	<i>height h, mm</i>	<i>neck</i>	<i>pack size</i>
clear glass					
23 180 17	50	44	78	24/20	10
23 180 24	100	54	100	29/22	10
23 180 36	250	72	130	34/24	10
23 180 44	500	89	165	45/27	10
23 180 54	1000	110	200	60/46	10
amber glass					
23 189 17	50	44	78	24/20	10
23 189 24	100	54	100	29/22	10
23 189 36	250	72	130	34/24	10
23 189 44	500	89	165	45/27	10
23 189 54	1000	110	200	60/46	10



Bottle, reagent, narrow neck, standard ground joint, DURAN borosilicate

<i>Cat. No</i>	<i>capacity ml</i>	<i>dia. d, mm</i>	<i>height h, mm</i>	<i>neck</i>	<i>pack size</i>
clear glass					
21 164 14	25	36	64	12/21	10
21 164 17	50	42	80	14/15	10
21 164 24	100	52	96	14/15	10
21 164 36	250	70	133	19/26	10
21 164 44	500	86	164	24/29	10
21 164 54	1000	107	200	29/32	10
21 164 63	2000	134	248	29/32	10
21 164 73	5000	182	324	45/40	1
21 164 86	10000	227	398	60/46	1
21 164 91	20000	288	492	60/46	1
amber glass					
21 167 14	25	36	64	12/21	10
21 167 17	50	42	80	14/15	10
21 167 24	100	52	96	14/15	10
21 167 36	250	70	133	19/26	10
21 167 44	500	86	164	24/29	10
21 167 54	1000	107	200	29/32	10
21 167 63	2000	134	248	29/32	10
21 167 73	5000	182	324	45/40	1
21 167 86	10000	227	398	60/46	1
21 167 91	20000	288	492	60/46	1



Bottle, reagent, narrow neck, with glass flat-head stopper, DURAN borosilicate


<i>Cat. No</i>	<i>capacity ml</i>	<i>dia. d, mm</i>	<i>height h, mm</i>	<i>neck</i>	<i>pack size</i>
clear glass					
21 165 08	10	28	52	10/19	10
21 165 14	25	36	64	12/21	10
21 165 17	50	42	80	14/15	10
21 165 24	100	52	96	14/15	10
21 165 36	250	70	133	19/26	10
21 165 44	500	86	164	24/29	10
21 165 54	1000	107	200	29/32	10
21 165 63	2000	134	248	29/32	10

SCHOTT LABORATORY GLASSWARE PRODUCT LISTING

Bottle, reagent, narrow neck, with glass flat-head stopper, DURAN borosilicate

<i>Cat. No</i>	<i>capacity ml</i>	<i>dia. d, mm</i>	<i>height h, mm</i>	<i>neck</i>	<i>pack size</i>
clear glass, continued					
21 165 73	5000	182	324	45/40	1
21 165 86	10000	227	398	60/46	1
21 165 91	20000	288	492	60/46	1
amber glass					
21 168 08*	10	28	52	10/19	10
21 168 14	25	36	64	12/21	10
21 168 17	50	42	80	14/15	10
21 168 24	100	52	96	14/15	10
21 168 36	250	70	133	19/26	10
21 168 44	500	86	164	24/29	10
21 168 54	1000	107	200	29/32	10
21 168 63	2000	134	248	29/32	10
21 168 68*	3000				1
21 168 73	5000	182	324	45/40	1
21 168 86	10000	227	398	60/46	1
21 168 91	20000	288	492	60/46	1


*until stocks sold



Bottle, reagent, narrow neck, with polyethylene stopper, DURAN borosilicate

<i>Cat. No</i>	<i>capacity ml</i>	<i>dia. d, mm</i>	<i>height h, mm</i>	<i>neck</i>	<i>pack size</i>
clear glass					
21 160 08	10	28	52	10/19	10
21 160 14	25	36	64	12/21	10
21 160 17	50	42	80	14/15	10
21 160 24	100	52	96	14/15	10
21 160 36	250	70	133	19/26	10
21 160 44	500	86	164	24/29	10
21 160 54	1000	107	200	29/32	10
21 160 63	2000	134	248	29/32	10
21 160 73	5000	182	324	45/40	1
21 160 86	10000	227	398	60/46	1
21 160 91	20000	288	492	60/46	1
amber glass					
21 169 08	10	28	52	10/19	10
21 169 14	25	36	64	12/21	10
21 169 17	50	42	80	14/15	10
21 169 24*	100	52	96	14/15	10
21 169 36	250	70	133	19/26	10
21 169 44	500	86	164	24/29	10
21 169 54	1000	107	200	29/32	10
21 169 63	2000	134	248	29/32	10
21 169 73	5000	182	324	45/40	1
21 169 86	10000	227	398	60/46	1
21 169 91	20000	288	492	60/46	1

* until sold



SCHOTT LABORATORY GLASSWARE PRODUCT LISTING


Bottle, reagent, narrow neck, with standard ground joint, SODA-LIME

<i>Cat. No</i>	<i>capacity ml</i>	<i>dia. d, mm</i>	<i>height h, mm</i>	<i>neck</i>	<i>pack size</i>
clear glass					
23 164 17	50	44	78	14/15	10
23 164 24	100	54	100	14/23	10
23 164 36	250	72	130	19/26	10
23 164 44	500	89	165	24/29	10
23 164 54	1000	110	200	29/32	10
23 164 63	2000	135	248	29/32	10
amber glass					
23 167 17	50	44	78	14/15	10
23 167 24	100	54	100	14/23	10
23 167 36	250	72	130	19/26	10
23 167 44	500	89	165	24/29	10
23 167 54	1000	110	200	29/32	10
23 167 63	2000	135	248	29/32	10




Bottle, reagent, narrow neck, with flat-head glass stopper, SODA-LIME

<i>Cat. No</i>	<i>capacity ml</i>	<i>dia. d, mm</i>	<i>height h, mm</i>	<i>neck</i>	<i>pack size</i>
clear glass					
23 165 17	50	44	78	14/15	10
23 165 24	100	54	100	14/23	10
23 165 36	250	72	130	19/26	10
23 165 44	500	89	165	24/29	10
23 165 54	1000	110	200	29/32	10
23 165 63	2000	135	248	29/32	10
amber glass					
23 168 17	50	44	78	14/15	10
23 168 24	100	54	100	14/23	10
23 168 36	250	72	130	19/26	10
23 168 44	500	89	165	24/29	10
23 168 54	1000	110	200	29/32	10



Bottle, reagent, narrow neck, with polyethylene stopper, SODA-LIME

<i>Cat. No</i>	<i>capacity ml</i>	<i>dia. d, mm</i>	<i>height h, mm</i>	<i>neck</i>	<i>pack size</i>
clear glass					
23 160 17*	50	44	78	14/15	10
23 160 24	100	54	100	14/23	10
23 160 36	250	72	130	19/26	10
23 160 44	500	89	165	24/29	10
23 160 54	1000	110	200	29/32	10
23 160 63*	2000	135	248	29/32	1
amber glass					
23 169 17	50	44	78	14/15	10
23 169 24	100	54	100	14/23	10
23 169 36	250	72	130	19/26	10
23 169 44	500	89	165	24/29	10
23 169 54	1000	110	200	29/32	10
23 169 63	2000	135	248	29/32	10




* whilst stocks last

SCHOTT LABORATORY GLASSWARE PRODUCT LISTING


Bottle, square, narrow neck, complete with dust cap and pouring ring, SODA-LIME

<i>Cat. No</i>	<i>capacity ml</i>	<i>DIN thread GL</i>	<i>diameter d, mm</i>	<i>height h, mm</i>	<i>pack size</i>
clear glass					
23 810 245	100	32 H	49	119	10
23 810 365	250	32 H	64	155	10
23 810 445	500	32 H	77	186	10
23 810 545	1000	45 H	96	223	10
amber glass					
23 816 245	100	32 H	49	119	10
23 816 365	250	32 H	64	155	10
23 816 445	500	32 H	77	186	10
23 816 545	1000	45 H	96	223	10




Bottle, square, wide neck, with DIN thread, without cap and pouring ring, SODA-LIME

<i>Cat. No</i>	<i>capacity ml</i>	<i>DIN thread GL</i>	<i>diameter d, mm</i>	<i>height h, mm</i>	<i>pack size</i>
clear glass					
23 820 24	100	32 K	49	109	10
23 820 36	250	45 K	64	146	10
23 820 44	500	54 K	77	173	10
23 820 54	1000	60 K	97	213	10
amber glass					
23 826 24	100	32 K	49	109	10
23 826 36	250	45 K	64	146	10
23 826 44	500	54 K	77	173	10
23 826 54	1000	60 K	97	213	10




Bottle, reagent, round, complete with cap and pouring ring, AMBER SODA-LIME

<i>Cat. No</i>	<i>capacity ml</i>	<i>DIN thread GL</i>	<i>diameter d, mm</i>	<i>height h, mm</i>	<i>pack size</i>
23 835 665	2500	45 H	138	285	1




Screw cap, tamper-proof, blue PP, for narrow neck square bottles

<i>Cat. No</i>	<i>DIN thread GL</i>	<i>diameter d, mm</i>	<i>height h, mm</i>	<i>pack size</i>
29 301 19	32 H	45	32	10
29 301 28	45 H	60	35	10



Screw cap, safety tamper-proof, red PP, for narrow neck square bottles

<i>Cat. No</i>	<i>DIN thread GL</i>	<i>diameter d, mm</i>	<i>height h, mm</i>	<i>pack size</i>
29 302 19	32 H	45	32	10
29 302 28	45 H	60	35	10



SCHOTT LABORATORY GLASSWARE PRODUCT LISTING


Screw cap, tamper-proof, blue PP, for wide neck square bottles

Cat. No	DIN thread GL	diameter d, mm	height h, mm	pack size
29 303 19	32 K	44	23	10
29 303 28	45 K	58	27	10
29 303 32	54 K	69	29	10
29 303 35	60 K	78	29	10



Dropping bottle, with interchangeable NS pipette complete with rubber teat, SODA-LIME

Cat. No.	capacity ml	height h, mm	pack size
clear glass			
23 270 17	50	79	10
23 270 24	100	105	10
amber glass			
23 270 176	50	79	10
23 270 246	100	105	10



Dropping bottle, with flat grooved stopper, clear, SODA-LIME


Cat. No.	capacity ml	pack size
23 272 15*	30	10

* whilst stocks last




Aspirator bottle, NS tubulated, without stoppers, DURAN borosilicate

Cat. No	capacity ml	NS joint	diameter d, mm	height h, mm	pack size
24 701 44	500	19/26	89	165	10
24 701 54	1000	19/26	110	200	10
24 701 63	2000	19/26	136	245	10
24 701 73	5000	29/32	186	318	1
24 701 86	10000	29/32	234	398	1
24 701 91	20000	29/32	300	492	1



Aspirator bottle, NS tubulated, complete with NS stopper and NS stopcock, DURAN borosilicate

Cat. No	capacity ml	NS joint neck	diameter d, mm	height h, mm	pack size
24 702 44	500	24/29	89	165	10
24 702 54	1000	29/32	110	200	10
24 702 63	2000	29/32	136	245	10
24 702 73	5000	45/40	186	318	1
24 702 86	10000	60/46	234	398	1
24 702 91	20000	60/46	300	492	1



Stopcock, NS, for aspirator bottles


Cat. No.	capacity ml	NS	pack size
24 148 03	500 - 2000	19/26	1
24 148 04	5000 - 20000	29 /32	1



SCHOTT LABORATORY GLASSWARE PRODUCT LISTING


Aspirator bottle, neck with GL45 screw thread, complete with GL32 tubulature, screw thread adapter coupling, silicone gasket and drain cock with ptfе spindle, DURAN borosilicate

Cat. No.	capacity ml	diameter d, mm	height h, mm	pack size
24 703 54	1000	101	225	1
24 703 63	2000	136	262	1
24 703 73	5000	183	335	1
24 703 86	10000	230	410	1




Draincock with ptfе spindle, for GL 32 screwthread for aspirator bottle

Cat. No.	bore mm	capacity ml	pack size
24 147 03	4	1000 + 2000	1
24 147 04	8	5000 + 10000	1




Desiccator, with flat flange and NOVUS NS junction. Outlet in lid with stopcock, DURAN borosilicate

Cat. No.	DN	flange o.d. mm	height h, mm	joint	pack size
24 782 46	100	153	174	24/29	1
24 782 57	150	215	239	24/29	1
24 782 61	200	270	296	24/29	1
24 782 66	250	320	344	24/29	1
24 782 69	300	380	420	24/29	1




Desiccator, with flat flange, without junction, with knobbed lid, DURAN borosilicate

Cat. No.	DN	flange o.d. mm	height h, mm	pack size
24 781 46	100	153	187	1
24 781 57	150	215	252	1
24 781 61	200	270	309	1
24 781 66	250	320	357	1
24 781 69	300	380	433	1




Desiccator base, with plane flange, DURAN borosilicate

Cat. No.	DN	flange o.d. mm	height h, mm	height h 1, mm	pack size
24 770 46	100	153	112	58	1
24 770 57	150	215	154	81	1
24 770 61	200	270	202	115	1
24 770 66	250	320	235	120	1
24 770 69	300	380	283	150	1



Desiccator base, with ring-grooved flange, DURAN borosilicate


Cat. No.	DN	flange o.d. mm	height h, mm	height h 1, mm	pack size
24 773 61	200	270	202	115	1



SCHOTT LABORATORY GLASSWARE PRODUCT LISTING


Desiccator base, with plane flange, with NS side outlet, NOVUS type, DURAN borosilicate

Cat. No.	DN	flange o.d. mm	height h, mm	height h ₁ , mm	pack size
24 771 46	100	153	112	58	1
24 771 57	150	215	154	81	1
24 771 61	200	270	202	115	1
24 771 66	250	320	235	120	1
24 771 69	300	380	283	150	1



Desiccator base, with plane flange, screw thread outlet, MOBILEX type, DURAN borosilicate

Cat. No.	DN	flange o.d. mm	height h, mm	height h ₁ , mm	pack size
24 772 57	150	215	154	81	1
24 772 61	200	270	202	115	1
24 772 66	250	320	235	120	1
24 772 69	300	380	283	150	1



Desiccator lid, with knob, DURAN borosilicate

Cat. No.	DN	flange o.d. mm	height h, mm	pack size
24 410 46	100	153	75	1
24 410 57	150	215	98	1
24 410 61	200	270	107	1
24 410 66	250	320	122	1
24 410 69	300	380	150	1



Desiccator lid, NOVUS type, with outlet for NS stopcock, DURAN borosilicate

Cat. No.	DN	flange o.d. mm	height h, mm	joint	pack size
24 420 46	100	153	62	24/29	1
24 420 57	150	215	85	24/29	1
24 420 61	200	270	94	24/29	1
24 420 66	250	320	109	24/29	1
24 420 69	300	380	137	24/29	1



Desiccator lid with special outlet, WERTEX, NS 24/29, DURAN borosilicate

Cat. No.	DN	flange o.d. mm	height h, mm	pack size
24 430 57	150	215	85	1
24 430 66	250	320	109	1



Desiccator lid with screw thread outlet, MOBILEX type, DURAN borosilicate

Cat. No.	DN	flange o.d. mm	height h, mm	GL	pack size
24 440 57	150	215	85	32	1
24 440 61	200	270	94	32	1
24 440 66	250	320	109	32	1
24 440 69	300	380	137	32	1



SCHOTT LABORATORY GLASSWARE PRODUCT LISTING

O-ring for 200mm desiccator base 24 773 61


Cat. No.	DN	d mm	s mm	pack size
----------	----	---------	---------	-----------

Perbunan (NBR), to withstand up to 100°C

29 214 61 200 236 5.3 1

Silicone rubber (VMQ), to withstand up to 230°C


29 215 61 200 236 5.3 1



Screw cap, GL32 red, with aperture for MOBILEX screw thread outlet

Cat. No.	d mm	h mm	pack size
----------	---------	---------	-----------

29 227 08 41 26 1



Silicone rubber sealing ring, GL32, for MOBILEX screw thread outlet

Cat. No.	pack size
----------	-----------

29 236 12 10



Screw cap, GL32, with seal, for MOBILEX screw thread junction

Cat. No.	d mm	h mm	pack size
----------	---------	---------	-----------

29 240 19 41 26 1



Desiccator plate

Cat. No	DN	d mm	pack size
---------	----	---------	-----------

porcelain

29 725 46 100 90 1

29 725 57 150 140 1

29 725 61 200 185 1

29 725 66 250 235 1

29 725 69 300 280 1

metal

29 080 46 100 90 1

29 080 57 150 140 1

29 080 61 200 185 1

29 080 66 250 235 1


29 080 69 300 280 1



Stopcock with ptfе spindle, for NOVUS side outlet

Cat. No.	NS	l mm	d mm	pack size
----------	----	---------	---------	-----------


24 798 03 24/29 108 8 1



Stopcock with ptfе spindle, for NOVUS tubulated lid

Cat. No.	NS	l mm	d mm	pack size
----------	----	---------	---------	-----------


24 799 04 24/29 124 8 1



SCHOTT LABORATORY GLASSWARE PRODUCT LISTING


Stopcock with ptfе spindle, for MOBILEX screw thread outlet

<i>Cat. No.</i>	<i>l</i> <i>mm</i>	<i>d</i> <i>mm</i>	<i>pack size</i>
24 797 03	160	8	1



Stopper, NS 24/29 for WERTEX safety outlet

<i>Cat. No.</i>	<i>l</i> <i>mm</i>	<i>h</i> <i>mm</i>	<i>pack size</i>
24 796 03	70	76	1



Test tubes, DURAN borosilicate


<i>Cat. No.</i>	<i>d x h</i> <i>mm</i>	<i>wall thick.</i> <i>mm</i>	<i>pack size</i>
-----------------	---------------------------	---------------------------------	------------------

with beaded rim

26 130 01	8 x 70	0.8 - 1.0	pk. 100
26 130 03	10 x 75	0.8 - 1.0	pk. 100
26 130 06	10 x 100	0.8 - 1.0	pk. 100
26 130 08	12 x 75	0.8 - 1.0	pk. 100
26 130 11	12 x 100	0.8 - 1.0	pk. 100
26 130 13	14 x 130	0.8 - 1.0	pk. 100
26 130 16	16 x 130	1.0 - 1.2	pk. 100
26 130 21	16 x 160	1.0 - 1.2	pk. 100
26 130 23	18 x 180	1.0 - 1.2	pk. 100
26 130 26	20 x 150	1.0 - 1.2	pk. 100
26 130 28	20 x 180	1.0 - 1.2	pk. 100
26 130 33	25 x 150	1.0 - 1.2	pk. 50
26 130 36	25 x 200	1.0 - 1.2	pk. 50
26 130 38	30 x 200	1.0 - 1.4	pk. 50

without rim

26 131 01	8 x 70	0.8 - 1.0	pk. 100
26 131 03	10 x 75	0.8 - 1.0	pk. 100
26 131 06	10 x 100	0.8 - 1.0	pk. 100
26 131 08	12 x 75	0.8 - 1.0	pk. 100
26 131 11	12 x 100	0.8 - 1.0	pk. 100
26 131 13	14 x 130	0.8 - 1.0	pk. 100
26 131 16	16 x 130	1.0 - 1.2	pk. 100
26 131 21	16 x 160	1.0 - 1.2	pk. 100
26 131 23	18 x 180	1.0 - 1.2	pk. 100
26 131 26	20 x 150	1.0 - 1.2	pk. 100
26 131 28	20 x 180	1.0 - 1.2	pk. 100
26 131 33	25 x 150	1.0 - 1.2	pk. 50
26 131 36	25 x 200	1.0 - 1.2	pk. 50
26 131 38	30 x 200	1.0 - 1.4	pk. 50




Test tube, SUPREMAX, with beaded rim, will withstand up to 700°C

<i>Cat. No.</i>	<i>d x h</i> <i>mm</i>	<i>wall thick.</i> <i>mm</i>	<i>pack size</i>
26 150 21	16 x 160	1.0 - 1.5	pk. 100

SCHOTT LABORATORY GLASSWARE PRODUCT LISTING


Test tubes, FIOLAX, with beaded rim

<i>Cat. No.</i>	<i>d x h</i> <i>mm</i>	<i>wall thick.</i> <i>mm</i>	<i>pack size</i>
26 110 01	8 x 70	0.4 - 0.5	pk. 100
26 110 03	10 x 75	0.4 - 0.5	pk. 100
26 110 06	10 x 100	0.4 - 0.5	pk. 100
26 110 08	12 x 75	0.4 - 0.5	pk. 100
26 110 11	12 x 100	0.4 - 0.5	pk. 100
26 110 13	14 x 130	0.4 - 0.5	pk. 100
26 110 16	16 x 130	0.5 - 0.6	pk. 100
26 110 21	16 x 160	0.5 - 0.6	pk. 100
26 110 23	18 x 180	0.5 - 0.6	pk. 100
26 110 26	20 x 150	0.5 - 0.6	pk. 100
26 110 28	20 x 180	0.5 - 0.6	pk. 100
26 110 33	25 x 150	0.6 - 0.7	pk. 50
26 110 36	25 x 200	0.6 - 0.7	pk. 50
26 110 38	30 x 200	0.7 - 0.8	pk. 50




Centrifuge tube, round bottom, DURAN borosilicate

<i>Cat. No.</i>	<i>d x h</i> <i>mm</i>	<i>capacity</i> <i>ml</i>	<i>pack size</i>
21 601 10	12 x 100	6	50
21 601 11	16 x 100	12	50
21 601 14	24 x 100	25	10
21 601 17	34 x 100	50	10
21 601 24	40 x 115	80	10
21 601 26	44 x 100	80	10
21 601 36	56 x 147	250	10




Centrifuge tube, conical, DURAN borosilicate

<i>Cat. No.</i>	<i>d x h</i> <i>mm</i>	<i>capacity</i> <i>ml</i>	<i>pack size</i>
24 263 09	16 x 100	12	50




Centrifuge tube, conical, DURAN borosilicate

<i>Cat. No.</i>	<i>d x h</i> <i>mm</i>	<i>capacity</i> <i>ml</i>	<i>pack size</i>
21 611 14	24 x 100	25	10
21 611 17	34 x 100	50	10



Culture tube, without rim, for use with Kapsenberg caps, DURAN borosilicate

<i>Cat. No.</i>	<i>d x h</i> <i>mm</i>	<i>wall thick.</i> <i>mm</i>	<i>matching cap</i>	<i>pack size</i>
26 132 21	16 x 160	1.0 - 1.2	29 010 09	pk:100
26 132 23	18 x 180	1.0 - 1.2	29 010 11	pk:100



SCHOTT LABORATORY GLASSWARE PRODUCT LISTING

Culture tube, with DIN thread, PBT screw cap, DURAN borosilicate

Cat. No.	d x h mm	DIN thread GL	pack size
26 135 115	12 x 100	14	pk:50
26 135 215	16 x 160	18	pk:50
26 135 235	18 x 180	18	pk:50




Culture tube, disposable, with DIN thread, PP screw cap with rubber gasket, SODA-LIME

Cat. No.	d x h mm	DIN thread GL	wall thick. mm	pack size
23 175 115	12 x 100	14	1	pk:100
23 175 145	16 x 100	18	1	pk:100
23 175 215	16 x 160	18	1	pk:100
23 175 235	18 x 180	18	1	pk:100



Culture tube, disposable, without rim, SODA-LIME

Cat. No.	h x d mm	wall thick. mm	pack size
23 172 018	75 x 9.75	0.60	pk:814
23 172 036	75 x 10.00	0.60	pk:766
23 172 058	75 x 11.75	0.80	pk:550
23 172 086	75 x 11.75	0.55	pk:550
23 172 096	75 x 12.25	0.55	pk:500
23 172 099	75 x 12.25	0.80	pk:500
23 172 078	100 x 11.75	0.80	pk:550
23 172 119	100 x 12.25	0.80	pk:500
23 172 148	100 x 15.50	0.80	pk:310
23 172 109	100 x 15.75	0.90	pk:310
23 172 128	120 x 12.25	0.80	pk:500
23 172 188	150 x 10.0	0.80	pk:766
23 172 198	150 x 15.5	0.80	pk:310
23 172 219	160 x 15.50	0.80	pk:310




Sample tube, NMR

Cat. No.	length mm	o.d. mm	i.d. mm	camber mm	pack size
Economic					
23 170 011	178	4.97 + .05	4.20 +.05	0.07	pk:50
Professional					
23 170 021	178	4.97±.025	4.20±.025	0.03	pk:25
Scientific					
23 170 031	178	4.97±.013	4.20±.025	0.013	pk:1



Spot plate, Feigl, Soda-lime


Cat. No.	l mm	b mm	h mm	pack size
23 671 52	130	100	14	10



SCHOTT LABORATORY GLASSWARE PRODUCT LISTING


Petri dish, DUROPLAN, DURAN borosilicate

Cat. No.	d 1 x h 1 mm	d 2 mm	h mm	pack size
21 755 41	60 x 20	54.5	22	10
21 755 43	80 x 20	74	22	10
21 755 46	100 x 15	94	17	10
21 755 48	100 x 20	94	22	10
21 755 51	120 x 20	114	22	10
21 755 53	150 x 30	143	32	10




Petri dish, SODA-LIME

Cat. No.	d x h mm	pack size
23 755 39	40 x 12	10
23 755 40	60 x 15	10
23 755 42	80 x 15	10
23 755 45	100 x 10	10
23 755 46	100 x 15	10
23 755 48	100 x 20	10
23 755 51	120 x 20	10
23 755 52	150 x 25	10
23 755 56	180 x 30	10
23 755 59	200 x 30	10
23 755 61	200 x 50	10




Petri dish, pressed, DURAN borosilicate

Cat. No.	d x h mm	pack size
without sections		
21 754 46	100 x 15	10
21 754 48	100 x 20	10
with two sections		
21 750 48	100 x 20	10
with three sections		
21 753 48	100 x 20	10
with four sections		
21 752 48	100 x 20	10




Culture flask, Fernbach type, bulbous shape, DURAN borosilicate

Cat. No.	capacity ml	d mm	d 1 mm	h mm	pack size
21 501 43	450	117	29	100	10



Culture flask, Fernbach type, conical shape, DURAN borosilicate


Cat. No.	capacity ml	d mm	d 1 mm	h mm	pack size
21 511 62	1800	200	45	158	2



SCHOTT LABORATORY GLASSWARE PRODUCT LISTING


Culture flask, Kolle, oval neck, DURAN borosilicate

Cat. No.	capacity ml	l x b mm	h mm	d 1 mm	pack size
21 521 41	400	140x200	39	39 x 60	10




Culture flask, Roux, DURAN borosilicate

Cat. No.	capacity ml	l x b mm	h mm	d 1 mm	pack size
21 541 58	1200	260x122	56	32	10




Culture flask, Roux, conical neck, eccentric, DURAN borosilicate

Cat. No.	capacity ml	l x b mm	h mm	pack size
21 571 58	1200	275x132	56	10




Penicillin flask, DURAN borosilicate

Cat. No.	capacity ml	l x b mm	h mm	d 1 mm	pack size
21 551 71	4000	370x20	90	50	1




Culture bottle, straight rim, for Kapsenberg cap 29 010 11, DURAN borosilicate

Cat. No.	capacity ml	d mm	d 1 mm	h mm	pack size
21 421 17	50	40	18	107	10
21 421 24	100	40	18	150	10
21 421 32	200	50	18	175	10



Erlenmeyer flask, straight rim, for Kapsenberg cap 29 010 11, DURAN borosilicate

Cat. No.	capacity ml	d mm	d 1 mm	h mm	pack size
21 491 24	100	64	18	120	10




Kapsenberg caps, aluminium

Cat. No.	for neck d 1 mm	pack size
29 010 09	16	10
29 010 11	18	10



Culture media bottle, with straight rim, for glass cap, DURAN borosilicate

Cat. No.	capacity ml	d mm	d 1 mm	h mm	pack size
21 431 39	300	71	31	169	10
21 431 44	500	83	46	204	10
21 431 54	1000	105	46	237	10



SCHOTT LABORATORY GLASSWARE PRODUCT LISTING


Glass cap, DURAN borosilicate

Cat. No.	for neck d 1 mm	pack size
21 441 18	31	10
21 441 29	46	10




Square bottle, Breed-Demeter, DURAN borosilicate

Cat. No.	capacity ml	d mm	b mm	h mm	pack size
21 481 31	180	28	48	148	10




Culture media bottle, with beaded rim, DURAN borosilicate

Cat. No.	capacity ml	d mm	d 1 mm	h mm	pack size
21 451 24	100	50	29	115	10
21 451 39	300	71	42	167	10
21 451 44	500	83	42	206	10
21 451 54	1000	105	46	235	10
21 451 66	2500	150	50	313	1
21 451 73	5000	185	54	390	1




Culture flask, Erlenmeyer shape, straight neck, for metal cap, DURAN borosilicate

Cat. No.	capacity ml	d mm	d 1 mm	h mm	pack size
21 771 24	100	64	38	114	10
21 771 32	200	79	38	138	10
21 771 36	250	85	38	149	10
21 771 39	300	87	38	161	10
21 771 44	500	105	38	183	10
21 771 54	1000	131	38	232	10
21 771 63	2000	166	38	305	10




Culture media bottle, straight neck, for metal cap, DURAN borosilicate

Cat. No.	capacity ml	d mm	d 1 mm	h mm	pack size
21 773 24	100	50	38	125	10
21 773 39	300	71	38	170	10
21 773 44	500	83	38	208	10
21 773 54	1000	105	38	240	10



Culture flask, Fernbach type, conical shape, straight neck for metal cap, DURAN borosilicate

Cat. No.	capacity ml	d mm	d 1 mm	h mm	pack size
21 774 62	1800	200	38	175	2



SCHOTT LABORATORY GLASSWARE PRODUCT LISTING

Metal caps, suitable for flasks 21 771, 21 774, bottles 21 773

Cat. No.	for neck d 1 mm	pack size
----------	--------------------	-----------

stainless steel

29 012 24 38 10

aluminium, anodised blue

29 013 24 10




Roller bottle for cell cultures, with DIN thread, blue PP cap and pouring ring, DURAN borosilicate

Cat. No.	DIN thread GL	d mm	h mm	pack size
----------	------------------	---------	---------	-----------

21 772 685 45 110 285 2

21 772 865 45 110 450 1



Bottle with rolled flange, without closure, DURAN borosilicate

Cat. No.	capacity ml	d mm	d 1 mm	h mm	pack size
----------	----------------	---------	-----------	---------	-----------


21 461 17 50 38 17 101 10

21 461 24 100 45 17 135 10

21 461 36 250 57 17 182 10

21 461 44 500 74 17 218 10

21 461 54 1000 95 17 265 10



Bottle with rolled flange, with clamp closure, DURAN borosilicate

Cat. No.	capacity ml	d mm	d 1 mm	pack size
----------	----------------	---------	-----------	-----------

21 465 17 50 38 17 10

21 465 24 100 45 17 10

21 465 36 250 57 17 10

21 465 44 500 74 17 10

21 465 54 1000 95 17 10



Staining jar, Coplin, for 10 microscope slides 76 x 26mm, SODA-LIME

Cat. No.	dia. a mm	h mm	pack size
----------	--------------	---------	-----------


23 319 00 45 100 10



Staining dish, Hellendahl, for 16 microscope slides 76 x 26mm, SODA-LIME

Cat.No.	l mm	b mm	h mm	pack size
---------	---------	---------	---------	-----------


23 314 00 90 40 90 10



SCHOTT LABORATORY GLASSWARE PRODUCT LISTING


Staining dish, Hellendahl, for 16 microscope slides 76 x 26mm, SODA-LIME

Cat.No.	<i>l</i> mm	<i>b</i> mm	<i>h</i> mm	pack size
23 315 00	60	60	100	10




Staining dish, Schieffedecker, for 10 microscope slides 76 x 26mm, SODA-LIME

Cat.No.	<i>l</i> mm	<i>b</i> mm	<i>h</i> mm	pack size
23 316 00	90	70	40	10



Staining tray, for 10 microscope slides 76 x 26mm or any width up to 52mm, DURAN borosilicate

Cat.No.	<i>l</i> mm	<i>b</i> mm	<i>h</i> mm	pack size
21 317 00	88	40	70	10




Stainless steel handle for staining tray

Cat. No.	pack size
29 075 00	10




Glass box with cover, for staining tray

Cat.No.	<i>l</i> mm	<i>b</i> mm	<i>h</i> mm	pack size
23 318 00	108	90	70	10



Volumetric flask, Class A, with PP stopper, DURAN borosilicate


Cat. No.	capacity ml	accuracy ±ml	<i>h</i> x <i>d</i> mm	stopper	pack size
21 678 07	5	0.025	70 x 24	7/16	10
21 678 08	10	0.025		7/16	10
21 678 12	20	0.04		10/19	10
21 678 14	25	0.04	110 x 40	10/19	10
21 678 17	50	0.06	140 x 49	12/21	10
21 678 24	100	0.10	170 x 61	12/21	10
21 678 25	100	0.10	170 x 61	14/23	10
21 678 32	200	0.15	210 x 76	14/23	10
21 678 36	250	0.15	220 x 81	14/23	10
21 678 44	500	0.25	260 x 101	19/26	10
21 678 54	1000	0.40	300 x 127	24/29	10
21 678 63	2000	0.60	370 x 161	29/32	10
21 678 73	5000	1.20	470 x 217	34/35	1



SCHOTT LABORATORY GLASSWARE PRODUCT LISTING


Mixing cylinder, graduated, hexagonal base, with PP stopper, CLASS B, DURAN borosilicate

<i>Cat. No.</i>	<i>capacity ml</i>	<i>grad. ml</i>	<i>h x d mm</i>	<i>stopper</i>	<i>pack size</i>
21 618 08	10	0.2	156 x 16	10/19	10
21 618 14	25	0.5	191 x 22	14/23	10
21 618 17	50	1.0	223 x 25	19/26	10
21 618 24	100	1.0	286 x 29	24/29	10
21 618 36	250	2.0	364 x 39	29/32	10
21 618 44	500	5.0	395 x 54	34/35	10
21 618 54	1000	10.0	500 x 66	45/40	1
21 618 63	2000	20.0	537 x 85	45/40	1




Measuring cylinder, graduated, hexagonal base, CLASS B, DURAN borosilicate

<i>Cat. No.</i>	<i>capacity ml</i>	<i>grad. ml</i>	<i>h x d mm</i>	<i>pack size</i>
21 396 07	5	0.1	112 x 13	10
21 396 08	10	0.2	137 x 14	10
21 396 14	25	0.5	167 x 21	10
21 396 17	50	1.0	196 x 25	10
21 396 24	100	1.0	256 x 29	10
21 396 36	250	2.0	331 x 39	10
21 396 44	500	5.0	360 x 53	10
21 396 54	1000	10.0	460 x 65	1
21 396 63	2000	20.0	500 x 85	1




Measuring cylinder, graduated, low form, CLASS B, DURAN borosilicate

<i>Cat. No.</i>	<i>capacity ml</i>	<i>grad. ml</i>	<i>h x d mm</i>	<i>pack size</i>
21 395 08	10	1.0	90 x 21	10
21 395 14	25	1.0	115 x 25	10
21 395 17	50	2.0	145 x 29	10
21 395 24	100	2.0	165 x 39	10
21 395 36	250	10.0	195 x 54	10
21 395 44	500	10.0	250 x 65	10
21 395 54	1000	20.0	285 x 85	1
21 395 63	2000	50.0	340 x 105	1



Burette with straight stopcock, glass key, CLASS B, DURAN borosilicate


<i>Cat. No.</i>	<i>capacity ml</i>	<i>grad. ml</i>	<i>h mm</i>	<i>pack size</i>
24 328 27	10	0.02	750	2
24 328 33	25	0.05	750	2
24 328 36	50	0.1	750	2
24 328 39	100	0.2	750	2



SCHOTT LABORATORY GLASSWARE PRODUCT LISTING


Burette with straight stopcock, glass key, Schellbach, CLASS AS, DURAN borosilicate

<i>Cat. No.</i>	<i>capacity ml</i>	<i>grad. ml</i>	<i>h mm</i>	<i>pack size</i>
24 329 27	10	0.02	750	2
24 329 33	25	0.05	750	2
24 329 36	50	0.1	750	2




Micro burette with straight stopcock, glass key, Schellbach, CLASS AS, DURAN borosilicate

<i>Cat. No.</i>	<i>capacity ml</i>	<i>grad. ml</i>	<i>h mm</i>	<i>pack size</i>
24 320 11	1	0.01	475	1
24 320 16	2	0.01	550	1
24 320 22	5	0.02	700	1




Automatic burette, Pellet pattern, with side stopcock, glass key, Schellbach, with 2 litre reservoir and blowing ball, CLASS AS, DURAN borosilicate

<i>Cat. No.</i>	<i>capacity ml</i>	<i>grad. ml</i>	<i>h mm</i>	<i>pack size</i>
24 318 275	10	0.02	930	1
24 318 335	25	0.05	930	1
24 318 365	50	0.1	930	1



Burette bottle, with 29/32 neck, DURAN borosilicate

<i>Cat. No.</i>	<i>capacity ml</i>	<i>h mm</i>	<i>d mm</i>	<i>pack size</i>
clear glass				
21 159 63	2000	200	160	1
amber glass				
21 159 636	2000	200	160	1



Volumetric pipette, CLASS B, AR glass

<i>Cat. No.</i>	<i>capacity ml</i>	<i>accuracy ±ml</i>	<i>length mm</i>	<i>pack size</i>
24 337 01	1	0.01	270	12
24 337 02	2	0.015	330	12
24 337 07	5	0.02	400	12
24 337 08	10	0.03	440	12
24 337 12	20	0.05	510	6
24 337 14	25	0.05	520	6
24 337 17	50	0.08	540	6
24 337 24	100	0.12	585	6



SCHOTT LABORATORY GLASSWARE PRODUCT LISTING

Volumetric pipette, CLASS AS, AR glass

Cat. No.	capacity ml	accuracy ±ml	length mm	pack size
24 338 01	1	0.007	270	12
24 338 02	2	0.01	330	12
24 338 07	5	0.015	400	12
24 338 08	10	0.02	440	12
24 338 12	20	0.03	510	6
24 338 14	25	0.03	520	6
24 338 17	50	0.05	540	6
24 338 24	100	0.08	585	6



Pipette filler, FLIP, for pipettes 0.1 to 100ml, red rubber

Cat. No.	pack size
29 270 10	1



Pipette filler, Safety, for pipettes up to 25ml, red rubber

Cat. No.	pack size
29 271 10	1




Pipette filler, Safety, for pipettes all pipette sizes, red rubber

Cat. No.	pack size
29 271 20	1

Funnel with short stem


Cat. No.	d mm	h mm	h ₁ mm	d ₁ mm	pack size
DURAN borosilicate					
21 351 23	35	60	35	6	10
21 351 28	45	80	45	6	10
21 351 33	55	95	55	8	10
21 351 38	70	120	70	8	10
21 351 41	80	140	80	10	10
21 351 46	100	180	100	10	10
21 351 51	120	210	120	16	10
21 351 57	150	250	150	16	10
21 351 59	180	290	150	20	1
21 351 61	200	325	175	26	1
21 351 66	250	360	175	30	1
21 351 69	300	409	175	30	1
SODA-LIME					
23 351 18	30	55	30	6	10
23 351 23	35	60	35	6	10
23 351 24	40	75	40	6	10
23 351 28	45	80	45	6	10
23 351 32	50	85	50	8	10
23 351 33	55	95	55	8	10



SCHOTT LABORATORY GLASSWARE PRODUCT LISTING


Funnel with short stem

<i>Cat. No.</i>	<i>d</i> <i>mm</i>	<i>h</i> <i>mm</i>	<i>h 1</i> <i>mm</i>	<i>d 1</i> <i>mm</i>	<i>pack size</i>
SODA-LIME, continued					
23 351 34	60	100	60	8	10
23 351 38	70	120	70	8	10
23 351 41	80	140	80	10	10
23 351 46	100	180	100	10	10




Powder funnel

<i>Cat. No.</i>	<i>d</i> <i>mm</i>	<i>h</i> <i>mm</i>	<i>h 1</i> <i>mm</i>	<i>d 1</i> <i>mm</i>	<i>pack size</i>
DURAN borosilicate					
21 354 33	55	60	30	22	10
21 354 38	70	70	30	22	10
21 354 41	80	80	30	24	10
21 354 46	100	100	30	26	10
21 354 51	120	105	30	35	10
21 354 55	160	140	30	35	1
21 354 61	200	170	30	40	1




Powder funnel

<i>Cat. No.</i>	<i>d</i> <i>mm</i>	<i>h</i> <i>mm</i>	<i>h 1</i> <i>mm</i>	<i>d 1</i> <i>mm</i>	<i>pack size</i>
SODA-LIME					
23 354 33	55	60	30	22	10
23 354 38	70	70	30	22	10
23 354 41	80	80	30	24	10
23 354 46	100	100	30	26	10




Funnel, long stem, DURAN borosilicate

<i>Cat. No.</i>	<i>d</i> <i>mm</i>	<i>h</i> <i>mm</i>	<i>h 1</i> <i>mm</i>	<i>d 1</i> <i>mm</i>	<i>pack size</i>
21 353 33	55	190	150	6	10
21 353 38	70	200	150	6	10
21 353 41	80	210	150	6	10



Funnel, long stem, SODA-LIME


<i>Cat. No.</i>	<i>d</i> <i>mm</i>	<i>h</i> <i>mm</i>	<i>h 1</i> <i>mm</i>	<i>d 1</i> <i>mm</i>	<i>pack size</i>
23 353 33	55	190	150	6	10
23 353 38	70	200	150	6	10
23 353 41	80	210	150	6	10
23 353 46	100	240	150	9	10



SCHOTT LABORATORY GLASSWARE PRODUCT LISTING


Analytical funnel, DURAN borosilicate

Cat. No.	d mm	h mm	h 1 mm	d 1 mm	pack size
21 331 37	65	200	150	9	10
21 331 41	80	210	150	9	10
21 331 48	110	265	180	9	10




Funnel, ribbed, DURAN borosilicate

Cat. No.	d mm	h mm	h 1 mm	d 1 mm	pack size
21 352 38	70	125	70	8	10
21 352 41	80	140	80	10	10
21 352 46	100	180	100	10	10
21 352 57	150	266	150	16	10
21 352 61	200	326	175	26	1




Filter funnel, conical, DURAN borosilicate

Cat. No.	porosity	capacity ml	d/d 1 mm	d 2/h mm	pack size
25 8540 3	3	25	55/25	8/100	1
25 8540 4	4	25	55/25	8/100	1



Filter funnel, DURAN borosilicate

Cat. No.	porosity	capacity ml	plate dia. mm	d 1/h mm	pack size
25 8520 X	1 - 3	50	35	10/130	1
25 8520 4	4	50	35	10/130	1
25 8520 5	5	50	35	10/130	1
25 8521 X	1 - 3	75	45	10/130	1
25 8521 4	4	75	45	10/130	1
25 8521 5	5	75	45	10/130	1
25 8522 X	1 - 3	125	60	10/140	1
25 8522 4	4	125	60	10/140	1
25 8522 5	5	125	60	10/140	1
25 8523 X	1 - 3	500	95	22/200	1
25 8523 4	4	500	95	22/200	1
25 8523 5	5	500	95	22/200	1
25 8524 X	1 - 3	1000	120	22/320	1
25 8524 4	4	1000	120	22/320	1
25 8524 5	5	1000	120	22/320	1
25 8526 1	1	4000	175	22/320	1
25 8526 X	2 - 3	4000	175	30/400	1
25 8526 4	4	4000	175	30/400	1




Note: the required porosity must be inserted in the catalogue number in place of X in the catalogue number

SCHOTT LABORATORY GLASSWARE PRODUCT LISTING

Buchner funnel, DURAN borosilicate


Cat.No.	capacity ml	plate dia. mm	d1 mm	h mm	pack size
21 341 22	70	48	10	130	1
21 341 28	125	60	10	140	1
21 341 34	220	73	18	180	1
21 341 44	500	95	22	240	1
21 341 54	1000	120	22	280	1



Filter crucible, DURAN borosilicate

Cat. No.	porosity	capacity ml	d mm	pack size
25 8510 X	2 - 3	8	24	10
28 8510 4	4	8	24	10
25 8511 X	1 - 3	15	28	10
25 8511 4	4	15	28	10
25 8512 X	1 - 3	30	36	10
25 8512 4	4	30	36	10
25 8512 5	5	30	36	10
25 8513 X	1 - 3	50	46	10
25 8513 4	4	50	46	10
25 8513 5	5	50	46	10

Note: the required porosity must be inserted in place of X in the catalogue number




Rubber sleeve, EPDM, for filter crucibles, sterilisable up to 180°C

Cat. No.	d mm	pack size
29 201 14	26	10
29 201 21	33	10
29 201 26	41	10
29 201 31	49	10




Adapter for filter funnels, DURAN borosilicate

Cat. No.	d mm	d1 mm	h mm	pack size
24 316 16	27	10	110	10
24 316 22	34	10	110	10
24 316 26	41	10	125	10
24 316 32	50	10	132	10



Filter head, with PP funnel and 2 FKM seals, DURAN borosilicate


Cat. No.	capacity ml	for disc Ø mm	thread Ø mm	pack size
24 720 24	30	24	28	1
24 720 50	250	50	54	1
24 720 90	1000	90	95	1



SCHOTT LABORATORY GLASSWARE PRODUCT LISTING

Filter head with screw thread, DURAN borosilicate

<i>Cat. No.</i>	<i>capacity ml</i>	<i>thread Ø mm</i>	<i>pack size</i>
24 721 24	30	28	1
24 721 35	250	54	1
24 721 90	1000	95	1



Filter disc with glass rim, DURAN borosilicate

<i>Cat. No.</i>	<i>porosity</i>	<i>disc Ø mm</i>	<i>pack size</i>
25 2024 X	1 - 3	24	1
25 2024 4	4	24	1
25 2050 X	1 - 3	50	1
25 2050 4	4	50	1
25 2090 X	1 - 3	90	1
25 2090 4	4	90	1



Note: the required porosity must be inserted in place of X in the catalogue number


FKM seal

<i>Cat. No.</i>	<i>for disc Ø mm</i>	<i>pack size</i>
29 220 24	24	10
29 220 50	50	10
29 220 90	90	10



PP funnel for filter head


<i>Cat. No.</i>	<i>thread size mm</i>	<i>d mm</i>	<i>pack size</i>
29 221 24	28	10	1
29 221 50	54	12	1
29 221 90	95	18	1



Pipeline filter, DURAN borosilicate

<i>Cat. No.</i>	<i>porosity</i>	<i>plate Ø mm</i>	<i>d 1 mm</i>	<i>pack size</i>
25 8550 X	1 - 3	30	14	1
25 8550 4	4	30	14	1
25 8551 X	1 - 3	60	16	1
25 8551 4	4	60	16	1
25 8552 X	1 - 3	90	16	1
25 8552 4	4	90	16	1

Note: the required porosity must be inserted in place of X in the catalogue number




SCHOTT LABORATORY GLASSWARE PRODUCT LISTING

Immersion filter, DURAN borosilicate

Cat. No.	porosity	plate Ø mm	d 1 mm	h mm	pack size
25 8556 X	1 - 3	35	10	210	1
25 8556 4	4	35	10	210	1


Note: the required porosity must be inserted in place of X in the catalogue number



Gas distribution tube, DURAN borosilicate


Cat. No.	porosity	d mm	d 1 mm	h mm	pack size
25 8560 X	0 - 2	6	22	250	5
25 8561 X	0 - 2	9	25	250	5
25 8562 1	1	11	34	250	5

Note: the required porosity must be inserted in place of X in the catalogue number




Gas washing bottle, screw cap system with Drechsel head, without filter disc, DURAN borosilicate

Cat. No.	capacity ml	DIN thread GL	pack size
24 713 00	500	45	1




Gas washing bottle, screw cap system, with Drechsel head, with filter disc, DURAN borosilicate

Cat. No.	capacity ml	DIN thread GL	porosity	pack size
25 7040 1	500	45	1	1



Gas washing bottle, with sintered disc, with standard ground joint, DURAN borosilicate

Cat. No.	capacity ml	porosity	joint size NS	h/h 1 mm	pack size
25 7010 1	100	1	34/35	250/180	1
25 7020 1	250	1	45/40	250/160	1



Gas washing bottle, with sealed-in sintered disc, 29/32 standard ground joint and cap, DURAN borosilicate


Cat. No.	capacity ml	porosity	d mm	h/h 1 mm	pack size
25 7030 1	350	1	60	250/180	1



Micro immersion filter, for reverse filtration, DURAN borosilicate

Cat. No.	porosity	d 1 mm	d mm	h mm	pack size
25 8571 X	1 - 3	10	6	100	10
25 8571 4	4	10	6	100	10

Note: the required porosity must be inserted in place of X in the catalogue number




SCHOTT LABORATORY GLASSWARE PRODUCT LISTING

Micro filter candle, without tube, DURAN borosilicate

Cat. No.	porosity	d mm	h mm	pack size
25 8572 X	1 - 3	13	22	10
25 8572 4	4	13	22	10


Note: the required porosity must be inserted in place of X in the catalogue number



Micro filter candle, with narrow tube, DURAN borosilicate

Cat. No.	porosity	cartridge d1/h1 mm	d mm	h mm	pack size
25 8573 X	1 - 3	13/22	8	200	10
25 8573 4	4	13/22	8	200	10


Note: the required porosity must be inserted in place of X in the catalogue number



Micro filter funnel, DURAN borosilicate

Cat. No.	porosity	capacity ml	d mm	d 1 mm	pack size
25 8574 X	1 - 3	0.8	10	6	10
25 8574 4	4	0.8	10	6	10
25 8575 X	1 - 3	2	10	6	10
25 8575 4	4	2	10	6	10

Note: the required porosity must be inserted in place of X in the catalogue number



Micro filter funnel, Pregl, DURAN borosilicate

Cat. No.	porosity	capacity ml	d mm	d 1 mm	pack size
25 8576 1	1	4	10	6	10



Filter tube, Allihn, DURAN borosilicate

Cat. No.	porosity	capacity ml	d mm	d 1 mm	pack size
25 8532 X	1 - 3	30	20	9	10
25 8532 4	4	30	20	9	10

Note: the required porosity must be inserted in place of X in the catalogue number



Filter disc, sintered glass, edges not fused, centred, surfaces untreated, DURAN borosilicate

Cat. No.	porosity	d mm	pack size
25 1510 X	0 - 3	10	10
25 1510 4	4	10	10
25 1520 X	0 - 3	20	10
25 1520 4	4	20	10
25 1525 X	0 - 3	25	10
25 1525 4	4	25	10
25 1530 X	0 - 3	30	10
25 1530 4	4	30	10
25 1540 X	0 - 3	40	10
25 1540 4	4	40	10

Note: the required porosity must be inserted in place of X in the catalogue number



SCHOTT LABORATORY GLASSWARE PRODUCT LISTING

Filter disc, sintered glass, edges not fused, centred, surfaces untreated, DURAN borosilicate, contd.


<i>Cat. No.</i>	<i>porosity</i>	<i>d</i> <i>mm</i>	<i>pack size</i>
25 1550 0	0	50	5
25 1550 X	1 - 3	50	5
25 1550 4	4	50	5
25 1556 0	0	60	1
25 1560 X	1-Mar	60	1
25 1560 4	4	60	1
25 1570 X	0 - 3	70	1
25 1570 4	4	70	1
25 1580 X	0 - 3	80	1
25 1580 4	4	80	1
25 1590 X	0 - 3	90	1
25 1590 4	4	90	1
25 1620 X	0 - 3	120	1
25 1620 4	4	120	1



Note: the required porosity must be inserted in place of X in the catalogue number


Filter flasks and bottles with side-arm socket, DURAN borosilicate

<i>Cat. No.</i>	<i>capacity</i> <i>ml</i>	<i>d</i> <i>mm</i>	<i>d 1</i> <i>mm</i>	<i>h</i> <i>mm</i>	<i>pack size</i>
flask shape					
21 183 36	250	85	35	155	10
21 183 44	500	105	34	185	10
21 183 54	1000	135	45	230	10
21 183 63	2000	166	60	255	1
bottle shape					
21 193 68	3000	170	60	295	1
21 193 73	5000	185	70	360	1
21 193 86	10000	240	70	420	1
21 193 88	15000	255	75	500	1
21 193 91	20000	290	75	535	1



Filter flasks and bottles with glass hose connection, DURAN borosilicate


<i>Cat. No.</i>	<i>capacity</i> <i>ml</i>	<i>d</i> <i>mm</i>	<i>d 1</i> <i>mm</i>	<i>h</i> <i>mm</i>	<i>pack size</i>
flask shape					
21 201 24	100	64	24	105	10
21 201 36	250	85	34	155	10
21 201 44	500	105	34	185	10
21 201 54	1000	135	45	230	10
21 201 63	2000	166	60	255	1
bottle shape					
21 191 68	3000	170	60	295	1
21 191 73	5000	185	70	360	1
21 191 86	10000	240	70	420	1
21 191 88	15000	255	75	500	1
21 191 91	20000	290	75	535	1



SCHOTT LABORATORY GLASSWARE PRODUCT LISTING


Filter flasks and bottles with PP plastic hose connection, DURAN borosilicate

<i>Cat. No.</i>	<i>capacity ml</i>	<i>d mm</i>	<i>d 1 mm</i>	<i>h mm</i>	<i>pack size</i>
flask shape					
21 204 245	100	65	24	105	10
21 204 365	250	85	35	155	10
21 204 445	500	105	34	185	10
21 204 545	1000	135	45	230	10
21 204 635	2000	166	60	255	1
bottle shape					
21 194 685	3000	170	60	295	1
21 194 735	5000	185	70	360	1
21 194 865	10000	240	70	420	1
21 194 885	15000	255	75	500	1
21 194 915	20000	290	75	535	1




Assembly set, with removable plastic hose connection, screw and seals

<i>Cat. No.</i>	<i>pack size</i>
29 258 54	10




Guko set, rubber gaskets, EPDM, conical, for filter flasks, sterilisable to 150°C 8 gaskets, from 21mm to 89mm

<i>Cat. No.</i>	<i>pack size</i>
29 202 00	1




Guko rubber gasket, EPDM, conical for filtering flask, sterilisable to 150°C

<i>Cat. No.</i>	<i>a mm</i>	<i>b mm</i>	<i>d mm</i>	<i>h mm</i>	<i>pack size</i>
29 202 12	21	12	2.5	18	1
29 202 17	27	17	3.0	20	1
29 202 23	33	21	3.0	24	1
29 202 27	41	27.5	4.0	27	1
29 202 32	53	33	5.0	34	1
29 202 36	68	48	5.5	35	1
29 202 39	78	58	6.0	35	1
29 202 43	89	66	6.5	40	1



Filtering device, complete, DURAN borosilicate

<i>Cat. No.</i>	<i>porosity</i>	<i>pack size</i>
25 7104 3	3	1




SCHOTT LABORATORY GLASSWARE PRODUCT LISTING

Filtration apparatus, Witt, DURAN borosilicate


<i>Cat. No.</i>	<i>porosity</i>	<i>capacity ml</i>	<i>pack size</i>
25 7230 X	1 - 3	20000	1
25 7230 4	4	20000	1
25 7230 5	5	20000	1

Note: the required porosity must be inserted in place of X in the catalogue number




Filter apparatus, Witt type, complete with interchangeable lid, with side plastic hose connection, DURAN borosilicate

<i>Cat. No.</i>	<i>DN</i>	<i>h 1 mm</i>	<i>h 2 mm</i>	<i>NS</i>	<i>pack size</i>
24 730 46	100	160	76	29/32	1
24 730 57	150	200	102	29/32	1
24 730 61	200	300	126	29/32	1



Filter apparatus, Witt type, complete with interchangeable lid, with side plastic hose connection, wide rough ground tubulature, DURAN borosilicate

<i>Cat. No.</i>	<i>DN</i>	<i>h 1 mm</i>	<i>h 2 mm</i>	<i>NS</i>	<i>pack size</i>
24 731 46	100	160	84	45/40	1
24 731 57	150	200	112	45/40	1




Lid for filter apparatus, Witt type, DURAN borosilicate

<i>Cat. No.</i>	<i>DN</i>	<i>NS</i>	<i>pack size</i>
24 398 46	100	29/32	1
24 398 57	150	29/32	1
24 398 61	200	29/32	1




Lid for filter apparatus, Witt type, with wide rough-ground tubulature, DURAN borosilicate

<i>Cat. No.</i>	<i>DN</i>	<i>NS</i>	<i>pack size</i>
24 450 46	100	45/40	1
24 450 57	150	45/40	1



Water jet pump, with non-return valve, DURAN borosilicate


<i>Cat. No.</i>	<i>water flow rate l/hour</i>	<i>min. water pressure bar</i>	<i>h mm</i>	<i>pack size</i>
24 362 99	300	1.2	275	1



SCHOTT LABORATORY GLASSWARE PRODUCT LISTING


Water jet pump, polypropylene, with non return valve, hose connection, adapters 1/2" and 3/4"

Cat. No.	connection for hose i.d. mm	water flow rate l/hour	min. water pressure bar	h mm	pack size
29 250 01	9 - 12	170	1	235	1




Reaction vessel, flat flange, with groove, for vacuum use, DURAN borosilicate

Cat. No.	capacity ml	DN	flange o.d. mm	d/h mm	pack size
24 390 24	100	60	100	70/85	1
24 390 36	250	60	100	70/125	1
24 390 44	500	100	138	106/120	1
24 390 54	1000	100	138	106/205	1
24 390 63	2000	100	138	140/270	1
24 390 71	4000	150	184	200/290	1
24 390 76	6000	150	184	215/320	1
24 390 86	10000	150	184	240/410	1




Reaction vessel, round bottom, flat flange with groove, for vacuum use, DURAN borosilicate

Cat. No.	capacity ml	DN	flange o.d. mm	d/h mm	pack size
24 395 63	2000	100	138	165/215	1
24 395 71	4000	100	138	206/265	1
24 395 76	6000	100	138	236/295	1
24 395 86	10000	100	138	280/340	1
24 395 91	20000	100	138	350/410	1




Beaker, flat flange with groove, for vacuum use, DURAN borosilicate

Cat. No.	capacity ml	DN	flange o.d. mm	d/h mm	pack size
24 394 54	1000	120	158	130/125	1
24 394 63	2000	120	158	130/200	1
24 394 68	3000	120	158	130/290	1
24 391 54	1000	150	184	154/120	1
24 391 63	2000	150	184	154/200	1
24 391 68	3000	150	184	154/265	1



Lid, flat flange, 4 standard ground necks, for vacuum use, DURAN borosilicate


Cat. No.	DN	flange o.d. mm	centre neck	side sockets	pack size
24 392 34	60	100	29/32	2x 19/26 angled 1x 14/23 angled	1
24 392 46	100	138	29/32	3x 29/32 angled	1
24 392 47	100	138	29/32	2x 29/32 angled 1x 14/23	1



SCHOTT LABORATORY GLASSWARE PRODUCT LISTING


Lid, flat flange, 4 standard ground necks, for vacuum use, DURAN borosilicate, continued

Cat. No.	DN	flange o.d. mm	centre neck	side sockets	pack size
24 392 51	120	158	29/32	parallel 2x 29/32 angled 1x 14/23	1
24 392 57	150	184	29/32	parallel 3x 29/32	1
24 392 58	150	184	29/32	parallel 3x 29/32	1
24 392 59	150	184	29/32	angled 2x 29/32 angled 1x 14/23	1
24 392 60	150	184	45/40	parallel 3x 29/32 angled	1




Lid, flat flange, 4 standard ground necks, for vacuum use, DURAN borosilicate

Cat. No.	DN	flange o.d. mm	centre neck	side sockets	pack size
24 396 46	100	138	29/32	2x 29/32 angled 1x 14/23 angled	1




Lid, flat flange, with centre neck only, for vacuum use, DURAN borosilicate

Cat. No.	DN	flange o.d. mm	centre neck	pack size
24 398 46	100	138	29/32	1
24 398 51	120	158	29/32	1
24 398 57	150	184	29/32	1
24 398 61	200	292	29/32	1
24 450 46	100	138	45/40	1
24 450 57	150	184	45/40	1



O-rings, FEP coated


Cat. No.	for flange DN	d x s mm	pack size
29 222 34	60	75 x 4	1
29 222 46	100	110 x 4	1
29 222 51	120	133 x 4	1
29 222 57	150	157 x 5	1
29 222 61	200	215 x 5	1



SCHOTT LABORATORY GLASSWARE PRODUCT LISTING


O-ring, silicone rubber (VMQ)

<i>Cat. No.</i>	<i>for flange DN</i>	<i>d x s mm</i>	<i>pack size</i>
29 225 34	60	75 x 4	1
29 225 46	100	110 x 4	1
29 225 51	120	133 x 4	1
29 225 57	150	157 x 5	1
29 225 61	200	215 x 5	1




Quick release clamp, stainless steel, with 3 retaining clips

<i>Cat. No.</i>	<i>for flange DN</i>	<i>pack size</i>
29 071 34	60	1
29 071 46	100	1
29 071 51	120	1
29 071 57	150	1
29 071 61	200	1




Supporting device for reaction vessel, nichrome for rigid support of the lid or the reaction vessel, with 2 support rods

<i>Cat. No.</i>	<i>for flange DN</i>	<i>pack size</i>
29 073 46	100	1
29 073 57	150	1




Woulff bottle, with 3 necks, for vacuum use, DURAN borosilicate

<i>Cat.No.</i>	<i>capacity ml</i>	<i>neck NS</i>	<i>d mm</i>	<i>pack size</i>
24 709 44	500	19/26	87	1
24 709 54	1000	24/29	113	1
24 709 63	2000	29/32	135	1
24 709 73	5000	34/35	185	1



Woulff bottle, with 3 necks and bottom tubulature, for vacuum use, DURAN borosilicate


<i>Cat. No</i>	<i>capacity ml</i>	<i>neck NS</i>	<i>bottom NS</i>	<i>d mm</i>
24 710 44	500	19/26	19/26	87
24 710 54	1000	24/29	19/26	113
24 710 63	2000	29/32	19/26	135
24 710 73	5000	34/35	29/32	185



SCHOTT LABORATORY GLASSWARE PRODUCT LISTING


Dropping funnel, cylindrical, with glass key stopcock, graduated, DURAN borosilicate

<i>Cat. No.</i>	<i>capacity ml</i>	<i>d mm</i>	<i>h mm</i>	<i>socket NS</i>	<i>pack size</i>
24 122 17	50 x 1	9	310	19/26	1
24 122 24	100 x 2	9	330	19/26	1
24 122 36	250 x 5	10	415	29/32	1
24 122 44	500 x 10	10	480	29/32	1
24 122 54	1000 x 20	13	560	29/32	1




Dropping funnel, cylindrical, with standard ground cone, graduated, DURAN borosilicate

<i>Cat. No.</i>	<i>capacity ml</i>	<i>cone NS</i>	<i>socket NS</i>	<i>h mm</i>	<i>pack size</i>
24 124 20	50 x 1	14/23	19/26	220	1
24 124 25	100 x 2	14/23	19/26	240	1
24 124 24	100 x 2	19/26	19/26	240	1
24 124 36	250 x 5	24/29	29/32	320	1
24 124 37	250 x 5	29/32	29/32	320	1
24 124 44	500 x 10	24/29	29/32	400	1
24 124 46	500 x 10	29/32	29/32	400	1
21 124 56	1000 x 20	29/32	29/32	480	1




Dropping funnel, pressure equalising, cylindrical, graduated, with standard ground cone, glass key stopcock, DURAN borosilicate

<i>Cat. No.</i>	<i>capacity ml</i>	<i>cone NS</i>	<i>socket NS</i>	<i>h mm</i>	<i>pack size</i>
24 125 20	50 x 1	14/23	19/26	240	1
24 125 25	100 x 2	14/23	19/26	270	1
24 125 24	100 x 2	19/26	19/26	270	1
24 125 36	250 x 5	24/29	29/32	350	1
24 125 37	250 x 5	29/32	29/32	380	1
24 125 44	500 x 10	24/29	29/32	430	1
24 125 46	500 x 10	29/32	29/32	430	1



Separating funnel, with glass key stopcock and NS stopper, DURAN borosilicate


<i>Cat. No.</i>	<i>capacity ml</i>	<i>stopper NS</i>	<i>h mm</i>	<i>d mm</i>	<i>pack size</i>
24 291 36	250	24/20	235	90	1
24 291 44	500	24/29	275	114	1
24 291 54	1000	29/31	295	132	1
24 291 66	2500	45/40	370	182	1
24 291 73	5000	45/40	425	222	1
24 291 86	10000	45/40	490	286	1



SCHOTT LABORATORY GLASSWARE PRODUCT LISTING


Separating funnel, conical, with ptfе key stopcock and plastic stopper, *graduated*, DURAN borosilicate

Cat. No.	capacity ml	stopper NS	stem dia. mm	h mm	pack size
10 648 05	100	19/26	9	230	1
10 648 06	250	29/32	10	280	1
10 648 07	500	29/32	10	320	1
10 648 09	1000	29/32	13	380	1




Calcium chloride cylinder, DURAN borosilicate

Cat. No.	h x d mm	cone NS	pack size
21 570 42	240 x 40	29/32	1
21 570 48	315 x 54	34/35	1



Plastic cap, closed, red PBT, with PTFE-protected seal

Cat. No.	DIN thread GL	pack size
29 240 08	14	10
29 240 11	18	10
29 240 13	25	10
29 240 19	32	10
29 240 28	45	10




Seal, silicone rubber, ptfе-protected, for use with PBT caps

Cat. No.	DIN thread GL	pack size
29 248 08	14	10
29 248 11	18	10
29 248 13	25	10
29 248 19	32	10
29 248 28	45	10



Plastic cap with aperture, red PBT


Cat. No.	DIN thread GL	aperture d, mm	pack size
29 227 05	14	9.5	10
29 227 06	18	11	10
29 227 09	25	15	10
29 227 08	32	20	10
29 227 10	45	34	10



SCHOTT LABORATORY GLASSWARE PRODUCT LISTING


Silicone rubber sealing ring, with vulcanised-on ptfe washers, for caps with aperture

<i>Cat. No.</i>	<i>DIN thread GL</i>	<i>d x d 1 mm</i>	<i>for tube Ø mm</i>	<i>pack size</i>
29 234 06	14	12 x 6	5.5 - 6.5	10
29 235 06	18	16 x 6	5.5 - 6.5	10
29 235 08	18	16 x 8	7.5 - 9.0	10
29 235 10	18	16 x 10	9.0 - 11.0	10
29 237 08	25	22 x 8	7.5 - 9.0	10
29 237 10	25	22 x 10	9.0 - 11.0	10
29 237 12	25	22 x 12	11.0 - 13.0	10
29 236 10	32	29 x 10	9.0 - 11.0	10
29 236 12	32	29 x 12	11.0 - 13.0	10
29 236 14	32	29 x 14	13.0 - 15.0	10
29 236 16	32	29 x 16	15.0 - 17.0	10
29 236 18	32	29 x 18	17.0 - 19.0	10
29 238 26	45	42 x 26	25.0 - 27.0	10
29 238 32	45	42 x 32	31.0 x 33.0	10




Silicone rubber seal (septa) for piercing, for caps with aperture

<i>Cat. No.</i>	<i>DIN thread GL</i>	<i>Ø d mm</i>	<i>thickness mm</i>	<i>pack size</i>
29 246 05	14	12	2	100
29 246 06	18	16	2	100
29 246 09	25	22	2	100
29 246 08	32	29	2	100
29 246 10	45	42	3	100




Plastic hose connection, straight, PP, with silicone rubber seal, for GL14 screw cap with aperture

<i>Cat. No.</i>	<i>d mm</i>	<i>d 1 mm</i>	<i>pack size</i>
29 255 06	9	5	10



Plastic hose connection, bent, PP, with silicone rubber seal, for GL14 screw cap with aperture

<i>Cat. No.</i>	<i>d mm</i>	<i>d 1 mm</i>	<i>pack size</i>
29 247 05	8	4	10




Spare seal, silicone rubber, for plastic hose connection

<i>Cat. No.</i>	<i>pack size</i>
29 220 09	10

SCHOTT LABORATORY GLASSWARE PRODUCT LISTING

Plastic stopper, polyethylene

<i>Cat. No.</i>	<i>nominal size</i>	<i>pack size</i>
29 204 02	7/16	10
29 204 03	10/19	10
29 204 04	12/21	10
29 204 06	14/23	10
29 204 07	19/26	10
29 204 08	24/29	10
29 204 09	29/32	10
29 204 11	34/35	1
29 204 12	45/40	1
29 204 13	60/46	1
29 204 16	85/55	1



Plastic stopper, polyethylene, short form

<i>Cat. No.</i>	<i>nominal size</i>	<i>pack size</i>
29 203 06	14/15	10
29 203 08	24/20	10
29 203 09	29/22	10